

## Chromosome Conformation Capture (3C) (alternate technique) (PROT31)



**Alice Horton, Peter Fraser and Tom Sexton**  
Laboratory of Chromatin & Gene Expression  
The Babraham Institute  
Cambridge, UK

Email feedback to:  
[tom.sexton@bbsrc.ac.uk](mailto:tom.sexton@bbsrc.ac.uk)

### Introduction

This is an alternative protocol for 3C that has been adopted by the author's lab. Much of it is identical to the [previous version \(PROT5\)](#). The major differences are in the amounts of DNA used at different steps. We reliably get the same results as previously, but also a greater yield of 3C material. The original yields are ample for real-time PCR analysis, but greater yields are required if the 3C material is going to be processed further.

The 3C (Chromosome Conformation Capture) technique generates a population average measurement of juxtaposition frequency between any two genomic loci, thus providing information on their relative proximity in the nucleus (Dekker *et al.*, 2002). Cells are fixed with formaldehyde which forms DNA-protein and protein-protein cross-links between regions of the genome in proximity (see [figure 1](#)). Subsequent restriction enzyme digestion and intra-molecular ligation produces novel junctions between restriction fragments in proximity in the nucleus. Novel ligation products can be detected by PCR. We adapted the 3C assay (Dekker *et al.*, 2002) to determine the conformation of mouse chromosome 7 and in particular the co-localization of actively transcribed genes in transcription factories (Osborne *et al.*, 2004). The 3C assay can also be used to reveal proximity between active genes and distal genomic elements (Tolhuis *et al.*, 2002).

An important consideration in the interpretation of 3C data is the understanding that not all pairs of restriction fragments that provide a positive result (i.e. generate a novel PCR product) are necessarily engaged in a functional interaction in the nucleus. For example, compare results of Tolhuis *et al.*, 2002, with Carter *et al.*, 2002 in which the higher order structure of the mouse hbb locus was assayed by two different methods, 3C and RNA TRAP. Clearly, distal fragments can be cross-linked by formaldehyde simply because they are near each other in the nucleus, and presumably can "bump into" each other during the fixation process (Osborne *et al.*, 2004). Therefore fixation conditions are critical in the 3C assay since increased fixation leads to greater cross-linking resulting in the detection of chromatin fragments that may be in proximity in the nucleus but not necessarily engaged in a specific intermolecular interaction with implied function.

### Procedure

#### Fixation and digestion

1. Dissect anaemic mouse spleen and kidney and strain through 70 $\mu$ m strainer into a chilled petrie dish in cold D-MEM + 10% FBS. Transfer to 50ml Falcon and make up to 50ml with cold medium;
2. Centrifuge at 1300 rpm for 8 minutes at 4°C. Remove supernatant and resuspend cells in trace of liquid;
3. Make up to 40ml with room temperature medium + 2% formaldehyde (2.16ml 37% formaldehyde + 37.84ml medium), and fix for 5 minutes at room temperature on rocker;
4. Quench with 5.7ml cold 1M glycine, and centrifuge at 1300 rpm for 8 minutes at 4°C;
5. Wash with 50ml cold PBS and centrifuge at 1300 rpm for 8 minutes at 4°C;
6. Make up to 50ml in cold [lysis buffer](#) and incubate on ice for 30 minutes with occasional mixing;
7. Centrifuge at 1800 rpm for 5 minutes at 4°C and resuspend nuclei in trace of liquid;
8. Dilute two to three-fold with 1.2x NEB3, and take 1x10<sup>7</sup> nuclei aliquots. Pulse centrifuge to visually check nuclei amount, and raise each 1x10<sup>7</sup> nuclei aliquot in 500 $\mu$ l 1.2x NEB3;
9. Add 7.5 $\mu$ l 20% SDS and incubate for 1 hour at 37°C, 950 rpm;

10. Add 50µl 20% Triton-X100 and incubate for 1 hour at 37°C, 950 rpm;
11. Add 30µl/1500U Bgl II and incubate overnight at 37°C, 950 rpm. Optional: Add 10µl/500U Bgl II in the morning, and incubate for a further 1-2 hours at 37°C, 950 rpm.

## Ligation

12. Add 40µl 20% SDS and incubate for 25 minutes at 65°C, 950 rpm;
13. Transfer to 15ml Falcon with 7ml 1.1x [ligation buffer](#);
14. Add 375µl 20% Triton-X100 and incubate for 1 hour at 37°C, mixing occasionally;
15. Add 2µl/800U T4 DNA ligase and incubate for 4 hours in a 16°C water bath, then for 30 minutes at room temperature.  
Optional: An aliquot can be used as a no ligase control (all subsequent treatments are the same);
16. Add 90µl/900µg proteinase K and incubate overnight at 65°C.

## DNA purification

17. Cool to room temperature and add 3µl/300µg RNase A. Incubate at 37°C for 1 hour;
18. Transfer to 50ml Falcon, add 10ml phenol and vortex. Centrifuge at 3500 rpm for 15 minutes at room temperature;
19. Take upper phase and extract as previously with 7ml chloroform;
20. Take upper phase and add 700µl 2M NaOAc, pH 5.2, and 17.5ml ethanol. Precipitate for 1-3 hours at -20°C;
21. Centrifuge at 3500 rpm for 30 minutes at 4°C and remove supernatant;
22. Vortex pellet with 20ml 70% ethanol and centrifuge at 3500 rpm for 30 minutes at 4°C;
23. Air dry pellet for 5 minutes at room temperature and 3 minutes at 37°C (with lid loosely on), then add 100µl water and incubate for 2 hours at 65°C;
24. Rigorously pipette and fully resuspend DNA overnight at 37°C.

## 3C validation


25. Quantitate DNA by PicoGreen assay.

Check for 3C products by PCR (two rounds of 35 cycles with nested primers).

## Materials & Reagents

<b>lysis buffer</b>	10mM Tris-HCl, pH 8 (500µl 1M) 10mM NaCl (100 µl 5M) 0.2% NP-40 (100 µl) 1 tablet complete protease inhibitors
<b>ligation buffer</b>	33 mM Tris-HCl, pH 8 (1.65ml 1M per 50ml) 11mM MgCl <sub>2</sub> (550µl 1M per 50ml) 11mM DTT (550 µl 1M per 50ml) 1.1 mM ATP (1.1ml 50mM per 50ml)

## Figures



**Figure 1.**

An overview of the 3C technique. Fixation with formaldehyde is followed by BglII digestion and intra-molecular ligation. Cross-links are reversed and novel ligation products are detected by PCR (adapted from Figure 1A. Dekker *et al.*, 2002).

## References

1. Dekker J, Rippe K, Dekker M, Kleckner N (2002) Capturing Chromosome Conformation. *Science* **295**:1306-11.

2. Osborne CS, Chakalova L, Brown KE, Carter D, Horton A, Debrand E, Goyenechea B, Mitchell J, Lopes S, Reik W, Fraser P (2004) Active genes dynamically colocalize to shared sites of ongoing transcription. *Nature Genetics* **36**:1065-1071.
3. Tolhuis B, Palstra RJ, Splinter E, Grosveld F, de Laat W (2002) Looping and Interaction between Hypersensitive Sites in the Active  $\beta$ -globin Locus. *Molecular Cell* **10**:1453-1465.
4. Carter D, Chakalova L, Osborne CS, Dai YF and Fraser P. 2002. Long-range chromatin regulatory interactions in vivo. *Nature Genetics* **32**: 623-6.
5. Spilianakis, C. G., and Flavell, R. A. (2004). Long-range intrachromosomal interactions in the T helper type 2 cytokine locus. *Nat Immunol*, **5**:1017-1027.
6. Vakoc, C., Letting, D. L., Gheldof, N., Sawado, T., Bender, M. A., Groudine, M., Weiss, M. J., Dekker, J., and Blobel, G. A. (2005). Proximity among Distant Regulatory Elements at the beta-Globin Locus Requires GATA-1 and FOG-1. *Mol Cell*, **17**:453-462.

Comments page 0 of 0:

[Add a Comment](#)

There are currently **0** comments to display.

